

# Sectra PACS Plugin



October 30, 2017

Software platform v2.1 R5964

MEDVISO AB  
<http://www.medviso.com>  
Griffelvägen 3  
SE-224 67 Lund  
Sweden  
Tel: +46-76-183 6442

# Contents

<b>1</b>	<b>Background</b>	<b>1</b>
1.1	Conditions for use . . . . .	1
<b>2</b>	<b>Conventions and Abbreviations</b>	<b>3</b>
2.1	Typographic conventions . . . . .	3
<b>3</b>	<b>Installation overview</b>	<b>5</b>
3.1	System Requirements . . . . .	5
3.2	Overview . . . . .	5
3.3	Installing Segment . . . . .	6
3.4	Installing Microsoft Visual C++ 2010 Redistributable Package	6
3.5	Enable communication from Sectra PACS to Segment . . . . .	6
3.6	Configure the Sectra PACS workstation . . . . .	7
3.7	Configuring communication . . . . .	9
3.7.1	Set up DICOM node in Sectra PACS . . . . .	9
3.7.2	Set up connection from Segment . . . . .	9
<b>4</b>	<b>Using the Sectra PACS Plugin</b>	<b>11</b>
4.1	Send images from Sectra PACS to Segment . . . . .	11
4.2	Send images from Segment to Sectra PACS . . . . .	11



# 1 Background

The Sectra PACS Plugin is as an addition to the software Segment, Segment CMR or Segment CT that allows for direct access to images in Sectra PACS. In this manual the software Segment, Segment CMR or Segment CT are used interchangeable since they work similarly with the Sectra PACS Plugin. At all places this manual specifies Segment this can also be read as Segment CMR or Segment CT.

Using the Sectra PACS Plugin it is possible to with a simple right click to send images and open them in Segment for analysis. From Segment it is possible with a single click to send the images back to the Sectra PACS including measurements and analysis etc.

The plugin is sold and distributed by Medviso AB. For more details about the software Segment, please see Segment User Manual. For more details about Segment CMR, please see Segment CMR User Manual. For more information about Sectra PACS, please see documentation supplied by SECTRA Imtec AB (<http://www.sectra.se>).

## 1.1 Conditions for use


- **Caution:** Federal law restricts this device to sale by or on the order of professionals trained in medical image analysis.
- The software is a tool to direct access images in Sectra PACS from Segment software packages.



# 2 Conventions and Abbreviations

The following typographic conventions are used in this manual.

## 2.1 Typographic conventions

A	Key A at the keyboard.
Ctrl-A	Control key. Hold down Ctrl key and A simultaneously.
	Icon in toolbar.
*.mat	Filename extension.
C:/Program	Folder.
File	Menu, e.g. File menu.
File→Save As	Sub menu, e.g. under the File menu the item Save As is found.
<input type="button" value="Close"/>	Push/Toggle button in the graphical user interface.
<input type="radio"/> Endocardium	Radiobutton in the graphical user interface.
<input checked="" type="checkbox"/> Single frame	Checkbox in the graphical user interface.





# 3 Installation overview

Before attempting the installation, please see system requirements below.

## 3.1 System Requirements

In order to run the SECTRA PASCs plugin you need to have:

- Access to Sectra PACS IDS7, including log-in.
- A powerful enough computer to run Segment. For further details, please see Chapter 'System Requirements' in the Segment User Manual.
- For computers equipped with Windows XP you need to have Service Pack 3 installed.
- Supported operating systems are Windows XP 32 bit, Windows XP 64 bit, and Windows 7 64 bit.

Please note that the plugin has not yet been tested on Sectra PACS IDS5. If you are commercially interested in an IDS 5 version, please contact [sales@medviso.com](mailto:sales@medviso.com).

## 3.2 Overview

The installation consist of three major steps:

1. Install Segment on each client computer you want to run the Sectra PACS Plugin.
2. Enable communication between Segment and Sectra PACS.
3. Configure the Sectra PACS workstation to allow plugins.
4. Configure Sectra PACS and Segment to be able to send images back from Segment to Sectra PACS.

Please note that the above process needs to be repeated on each client computer that will run Segment and the Sectra PACS Plugin.

### 3.3 Installing Segment

Please follow the instructions on how to install Segment from Segment Installation Manual. If you want to deploy Segment to many machines, it is technically possible to perform a network installation of the Segment and the MCRInstaller described in the first step. This may however make Segment to start slower.

You will need to copy a license file to the folder where you have installed Segment. This license file is supplied to you either at the time of purchase of Sectra PACS Plugin, or by [sales@medviso.com](mailto:sales@medviso.com) for evaluation purposes. The license file is necessary to run the Sectra PACS Plugin.

### 3.4 Installing Microsoft Visual C++ 2010 Redistributable Package

To successfully register the plugin, currently the runtime components for Visual C++ libraries have to be installed. These can be downloaded from

<http://www.microsoft.com/download/en/details.aspx?id=5555> (32 bit)

<http://www.microsoft.com/download/en/details.aspx?id=14632> (64 bit)

### 3.5 Enable communication from Sectra PACS to Segment

This step will enable the communication from Sectra PACS to Segment.

This step will register an .dll file. If you are using 32 bit platform, then double click on the file `sectraplugin32.bat` in Segment installation folder. If you are running a 64 bit system, then click on `sectraplugin64.bat`. You will now see a message box that looks like in Figure 1.

If the installation was not successful please contact [support@medviso.com](mailto:support@medviso.com) for help with trouble shooting. Please note that you need to have Service Pack 3 installed on Windows XP.




Figure 1: Successful notification .dll installation.

## 3.6 Configure the Sectra PACS workstation

Before you can use the Sectra PACS Plugin you need to configure the PACS installation on each computer where you want to run the application. Depending on the PACS master settings this step may not be required. In order to perform this step you need to have PACS administration system privileges.

1. Start Sectra PACS.
2. Login in as local PACS administrator.
3. In the menu Tools or Verktyg select Options or Alternativ. The keyboard short cut for this is **Shift-Ctrl-0**. Select the name of the workstation as the role. Click continue.
4. Click on the tab for Advanced Settings.
5. Select CAI Settings.
6. Adjust settings according to Figure 2.
7. You are now set to go. However, we strongly recommend that you first test it, and then also logoff and test it for users with standard privileges.

## CHAPTER 3. INSTALLATION OVERVIEW

---




Figure 2: Local PACS settings for each client computer.

## 3.7 Configuring communication

This step is optional and only required if you want to be able to send images and reports to back Sectra PACS from Segment.

### 3.7.1 Set up DICOM node in Sectra PACS

Ask your local Sectra PACS Administrator to set up a DICOM node (port) for the Segment Application. As PACS administrator you need to decide on the following:

- **Called AE.** This should be the AE Title for the DICOM node (or for the entire Sectra PACS).
- **Peer IP.** This is the IP address of the DICOM node (or the entire Sectra PACS).
- **Peer Port.** This is the port on which the DICOM node (or the entire Sectra PACS) is listening for incoming association requests.
- **Port.** This is the port suggested by Segment to be used in the communication.

You need to collect the following information about each client that you need to set up.

- Hostname
- IP number. [Need to check if this is necessary...]

### 3.7.2 Set up connection from Segment

You should create a new PACS connection in Segment and enter the information from above. This is done in the Preferences menu, and by clicking on **PACS Settings**. Thereafter, select **Add** under the **Define PACS connection** tab. The new settings should look something as shown in Figure 3.

For further details on how to configure PACS connections, please see Segment User Manual.




Figure 3: Example of communication settings for Sectra PACS.

# 4 Using the Sectra PACS Plugin


In this chapter it is assumed that you are well acquainted with using Sectra PACS. If you need assistance using the Sectra PACS, please consult you local hospital PACS support, or contact Sectra Imtec AB.

## 4.1 Send images from Sectra PACS to Segment

To open images in Segment, first open them in Sectra PACS. In the image window, simply right click on the desired image series. The menu showed in Figure 4. In the menu, click on Clinical Applications, or Kliniska Applicationer and select Segment.

Segment will automatically start and load the desired image stack. If Segment is already open, then it will load it into the current image stack collection. If you are loading images into another patient a warning message will appear.

## 4.2 Send images from Segment to Sectra PACS

After image analysis and quantification in Segment it is possible to send images back to Sectra PACS. This is done in Segment with the icon . It is also possible to send images from the Segment PACs Connection Module, or the Segment Patient Database.

The transmitted images are in Segment internal .mat format, but with a DICOM wrapper. This will enable that the images are visible as a thumbnail in Sectra PACS. When they are opened in Sectra PACS all you will see is an image preview. When opened in Segment through the Sectra PACS Plugin the full image and all measurements will be available. The internal representation used in Segment is stored in the private tag (FB37,4719).




Figure 4: Right click menu in Sectra PACS.