

Segment version 2.1 is released

## SEGMENT v2.1 (R5739) IS NOW RELEASED

The new version contains the new major features:

- \* Major rework of user interface
- \* Larger and clearer icons. All icons are moved to one location
- \* New row of icons that depend on active mode
- \* More functions that are available from icons
- \* Clearer report panel
- \* Flow graph is displayed in report panel
- \* Import contours are automatically aligned (import cine contours to LGE for example)
- \* Improved flow unwrapping
- \* Additional hot keys
- \* ...

As always we love feedback. Please send us your feedback on the new interface. Are there things that are missing? How can it be further improved?

## SEGMENT COMMUNITY

We are preparing to release Segment in source code version on GitHub and provide a platform for sharing user developed plugins. More information will be released in the autumn. Please contact Einar Heiberg if you want to be a pilot in this program.

## SEGMENT NEWSLETTER

You are receiving this newsletter since you are registered as a user of the medical image analysis software Segment. The software Segment is brought to you by Medviso AB in close collaboration with Lund Cardiac MR Group, Lund University, Sweden and Region of Scania, Skåne University Hospital. Segment is freely available for research purposes, provided that you properly refer to the software according to instructions on <http://medviso.com/research/how-to-refer/>

Best regards

Einar Heiberg, CTO and founder Medviso AB  
einar@medviso.com