




SEGMENT CT



Segment CT is a CE-marked vendor independent software for analysis of cardiovascular CT images.




Segment CT provides tools for fast, accurate and reproducible analysis of cardiovascular CT images. The software provides automatic segmentation algorithm for the LV.




Viewing

- Ease of use
- Multi panel
- View in 4 orthogonal projections
- Reconstruction tools from transversal to short-axis slices



Quantitative analysis

- LV function
- RV function
- Measurements in 3D
- Region of interest (ROI) analysis
- Signal intensity analysis
- Bullseye analysis



Global and Regional function

- EF, EDV, ESV and LVM
- Wall thickness
- Wall thickening
- Bullseye plot based on 17 segment AHA model

Loading data

- Advanced DICOM reader
- Integrated patient database
- Integrated PACS connection
- Plugin to Sectra PACS

Report function

- Include optional quantification results
- Images
- Reference values
- Comments by reviewing doctor
- Adjustable for your clinic

REGULATORY STATUS

SEGMENT CT bears the CE marking of conformity and is certified according to the ISO 13485 standard.

Medviso AB was founded in 2007 by Associate Professor Einar Heiberg as a spin-off from research at Lund University Sweden.



www.medviso.com
info@medviso.com
Phone: +46-76-1836442

Medviso AB | Griffelvägen 3 | SE- 224 67 Lund | Sweden

FREE EVALUATION
LICENSE AT
WWW.MEDVISO.COM